



Remember to Clean and Replace Your Equipment

MASK FRAME

Clean Weekly | Replace every 3 months

HEADGEAR

Clean Weekly | Replace every 6 months

TUBING

Clean Weekly | Replace every 3 months

HUMIDIFER CHAMBER *(Use only distilled water)*

Empty Daily | Clean Weekly | Replace every 6 months

DISPOSABLE FINE FILTERS

Change every 2 weeks


FULL FACE CUSHIONS

Wipe Daily | Replace every month

NASAL CUSHIONS AND PILLOWS

Wipe Daily | Replace every 2 weeks

Your supplies should be replaced when they show signs of wear. Cleaning should be done using a mild soap that is not antibacterial (i.e., Ivory, Palmolive, or baby shampoo/wash). You must meet your insurance carrier's compliance requirements to order. Coverage will vary based on insurance plan. Copays and deductibles will apply.



Visit us at [careprohs.com](https://www.careprohs.com)

CarePro
HOME MEDICAL 

Connect with people and answers you need.